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In order to ensure reliable operation of Huber Temperature control machines, it is recommended that regular maintenance is carried out. Unexpected and sudden failures as well as operating costs can be reduced by maintaining the machine.

The work to be carried out during servicing depends on the type of machine. Service can either be carried out by a Huber service technician or, depending on the service work to be done, by an in-house technician.

Your benefits of regular maintenance:

- 🗑️ extended service life
- 🗑️ Obtain optimum performance from your temperature control system
- 🗑️ Downtime minimized
- 🗑️ Reduced operating costs due to machine being in optimum condition
- 🗑️ Maintenance costs transparent and can be planned into maintenance budget

When is a regular maintenance to be recommended?

- applications with high downtime costs
- large temperature control units (e.g. Unistat & Unichiller -standing model)
- units, which have been in use for several years without any service
- units, which are in continuous operation (24/7)

Simple maintenance (especially for table-top units) may be carried out by an operating technician. For more information regarding the maintenance execution, refer to your operating manual.

Preventive maintenance could include just cleaning the condenser, the filter, and testing the safety devices. This can be carried out by an in-house technician. A service technician who has experience with Huber units could provide a professional maintenance, e.g. with data logging, checking cooling power and performance of the unit, testing electrical and mechanical components etc.

Depending on the prevailing regulations, the refrigeration machine has to be additionally checked for leaks. This usually depends on the quantity and type of refrigerant. Please find further information on page 8.

► **MAINTENANCE REMINDER – Message on controller**

On units with **Pilot ONE**-controller you will receive a reminder message which is displayed 4 weeks before the date when maintenance is due. This message does not affect the function or operation of the unit.

The maintenance reminder allows better management of the unit (e.g. reminder for calibration, preventive maintenance).

The standard interval is set to 12 months and can be changed in the menu between 1 month and 5 years.

► **F-GASES-REGULATION (EC) No. 517/2014**

F gases regulation (EC) No. 517/2014 of April 16, 2014, on fluorinated greenhouse gases, and repealing Regulation (EC) No. 842/2006.

The new Ordinance on Fluorinated Greenhouse Gases contributes to reducing the emissions of the industrial sector by 70 percent by the year 2030 when compared to 1990. By 2030, the new regulations are to reduce the emissions of fluorinated greenhouse gases (F gases) within the EU by 70 million tons of CO₂ equivalent to 35 million tons of CO₂ equivalent. Three main regulatory approaches are to reduce the emission of fluorinated greenhouse gases:

- Implementation of a gradual limitation (phase down) of the volume of hydrofluorocarbons (HFCs) available in the market. By 2030, the sales volumes are to be reduced to one-fifth of today's volumes.
- Issuing legal prohibitions on utilizing and placing F gases on the market if technically feasible and environmentally friendlier alternatives are available.
- Maintaining and supplementing the provisions on leakage tests, certification, labeling and disposal.

These regulations deal with all systems that contain fluorinated refrigerants. The substances dealt with in the Regulation (EC) No. 1005/2009 of the European Parliament and of the Council of May 16, 2009, on substances that deplete the ozone layer, are excluded (CFC/HCFC).

The directive regulates the reduction of the emission, utilization, recovery, and destruction of certain fluorinated greenhouse gases. It also regulates the identification and disposal of products and devices that contain these gases. Since 4 July 2007, responsible bodies must regularly check their stationary refrigeration system for leaks and have any leaks eliminated immediately.

Directive (EC) No. 303/2008 contains stipulations on the training and certification of companies and personnel that are permitted to execute the specified activities.

Obligations of the responsible body:


- Directive (EC) No. 842/2006 already imposed a number of obligations upon responsible bodies regarding certain fluorinated greenhouse gases. The new Ordinance on Fluorinated Greenhouse Gases upholds these to a large extent. Some duties are added while others are designed differently by this new ordinance. Please refer to the text of this ordinance for a complete overview of the individual responsible bodies.
- General obligation to reduce emissions.
- Only certified companies may maintain, repair or decommission refrigeration systems. The responsible bodies must verify that these companies are certified.
- ***As of 01.01.2017: Regular leak tests of stationary refrigeration systems by certified personnel. The required test interval is based on the refrigerant filling capacity and the type of refrigerant, converted to CO2 equivalent.***
- Responsibility of plant responsible bodies to recover F-Gases by certified personnel.
- Obligatory documentation requirements in the refrigeration system's operation manual, specifying type and volume of refrigerant used or recovered. The responsible authority must retain this documentation for at least 5 years and be able to produce it upon request of the responsible body.
- Temperature control units with natural refrigerants (NR) are exempt from this Regulation.
- Please refer to the data sheet or name plate of your temperature control unit for the quantity and type of refrigerant.

F-Gases Check required?

The following procedure explains how to check if unit has to be tested according (EC) No. 842/2006.

1. Type plate (on the back of the unit) or data sheet

e.g. refrigerant type: R507
 refrigerant quantity: 4,1 kg

2. Check prevailling national Regulation (e.g. EU Nr.517/2014)
 e.g. R507 = limit 2,5 kg < 4,1 kg
 The system has to be checked for leaks annually.

► ACTUAL LIMITS (EU Nr.517/2014)

Following limits (January 2017) are showing if a F-Gases check is required or not. If the system has an higher refrigerant quantity, an annual leak test will be required.

R507	→	from 2,50 kg	R23	→	from 0,68 kg	R452A	→	ab 4,67 kg
R134a	→	from 7,00 kg	R404a	→	from 2,55 kg			
R508b	→	from 0,75 kg	R14	→	from 1,35 kg			
Isceon 89	->	from 2,63 kg	R449A	→	ab 7,16 kg			

The stages on multistage units can be treated separately like two different cooling systems.

► reset MAINTENANCE REMINDER

For confirming the maintenance or in case maintenance is not required, the pictured procedure has to be followed.

After confirming the maintenance via “Acknowledge” a code will be displayed. This code has to be entered.